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Safety Guidelines for The Use of Syringe Needles

NO. 022

1. When dispensing liquids, avoid applying excessive pressure to the piston, which may cause the needle to pop out.
2. Keep needles away from fingers and other parts of the body. It is forbidden to transfer sharps directly by hand. Pallets can be used to transfer sharps.
3. Place uncapped needles in closed secondary containers during experiment operation and transport. Syringes with needles can be placed in plastic centrifuge tubes on a rack.
4. Do not bend or break contaminated needles to avoid injury and aerosols.
5. Syringes larger than 10 ml should not be used to handle pyrophoric chemicals. The piston of large syringes can easily be ejected accidentally, resulting in the splashing of burning, spontaneous liquids on the individual.
6. Disposed needles should not be capped. They should be put directly into the sharps box and not mixed with other medical waste. Sharp boxes need to be closed when they are 75 % full for subsequent processing.
7. Remove the needle by using the sharp box drip hole or the plier, it is prohibited to remove the needle by hand.

8. If the needle must be capped, use a tool or a suitable method (e.g., "one-handed scoop") to avoid being stuck.

"One-handed scoop" Procedure

① Place cap on hard flat surface. Slide syringe needle into cap.

② Scoop cap up with syringe needle so that cap is sitting on needle.

③ Using a hard flat surface, press cap onto needle until cap snaps into place.