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Safety Moment

How To Safely Use The Hydrothermal Autoclave Reactor (I)


A hydrothermal reactor, also known as a digestion vessel or high-pressure pot, is composed of a stainless steel jacket and a corrosion-resistant liner, usually made of polytetrafluoroethylene (PTFE). It typically has a volume of 25-500mL and can withstand certain temperatures and pressures. Common operating temperatures for hydrothermal bombs are between 180-220°C, with a design pressure limit of 1-3MPa.

Pre-use:

1. Attend training
2. Wear PPE
3. Inspect the exterior (stainless steel jacket and liner) for cracks, pitting, rust, creep, excessive damage, PTFE liner distortion, steel shell rupture, defects, and inadequate cleaning; do not use mismatched autoclave accessories.
4. Ensure pressure and temperature are below product limits. Avoid overloading.
5. Install hydrothermal reactor properly. Don't tighten off the bench. No tilting
6. Do not exceed 2/3 capacity with reaction materials
7. Use ovens with over-temperature protection or alarms for heating. Do not use muffle furnaces with poor low-temperature control for hydrothermal reactor experiments

From Zhaorong Wan and HSE

Minimum PPE to Enter Lab

Additional PPE Appropriate for the Task as Determined by Risk Assessment
Copyright © 2020 Small World Initiative, Inc. Artist: Sophia Gaynes.

The following systems are strictly prohibited for hydrothermal reactor experiments:

- Containing radioactive materials
- Containing explosive substances
- Containing chemically unstable substances that may decompose or at set temperatures
- Contaminated with used needles
- Containing perchloric acid
- Mixtures of nitric acid and organic compounds