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SAFETY ALERT

MAR 14 NO. 006

Incident of resuspended E. coli strain liquid splashing into eyes

Brief Description:

On January 22, 2025, in the afternoon, when a student used a pipette to resuspend E. coli in a centrifuge bottle on the laboratory bench, 2 drops of E. coli expression liquid splashed into the right eye from the blockage of the pipette tip. The student immediately used the eyewash station on the sink to rinse her eyes and then went to the clinic for treatment immediately when she found her eyes swollen. The doctor rinsed her eyes with plenty of water and prescribed levofloxacin eye drops.

Analysis of Factors affecting the accident:

- 1. When the student used a pipette to resuspend E. coli, to see clearly the suspension in the centrifuge bottle, she faced the mouth of the bottle, but she accidentally touched the pipette tip, causing droplets to splash into her eye.
- 2. Insufficient safety awareness. The student was not dressed properly. She did not choose appropriate face shielding based on risk factors when performing any experiments that may cause liquid splashing.

Learning Points:

- 1. Follow the laboratory dress code. Please wear a lab coat and safety glasses when entering the laboratory; if liquid splashing may occur during experimental operations, please choose appropriate PPE according to risk factors to protect your mouth, eyes and face to avoid injuries.
- 2. Follow microbiological standard operating procedures. Handling biological reagents or samples in a biological safety cabinet, and lower the sash to 160-250mm from the cabinet's surface.
- 3. Be aware of the risk of liquid splashing. Any experimental operation that may cause liquid splashing must be performed carefully and slowly. The volume of the solution in the unsealed container should not exceed 2/3 of the container to reduce the probability of liquid splashing.
- 4. It is recommended to choose reagent bottles made of transparent material so it is easy to observe the conditions inside.
- 5. Understanding laboratory emergency procedures. When liquid splashes into the eyes during the experiment, it is recommended to flush the eye with clean water immediately for at least 15min, report it to the laboratory safety officer in a timely manner. For details, please refer to "the Hong Kong University of Science and Technology (Guangzhou) Laboratory Safety Emergency Guide".

Big Data Warning:

After big data retrieval, there have been many similar accidents on social platforms like some victims have corneal damage, and in serious cases, even permanent vision loss.

Reference

NO. 017 Biosafety (II) - HSE HKUST GZ EN 第十三章 个人防护装备 - HSE HKUST GZ EN Emergency, Alert and Response Guidelines in Laboratories 洗眼设备和应急喷淋程序