SAFETY

ALERT

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JAN 23
NO. 005

Fire Incident Involving Laboratory “Oil Bath”

Incident Overview

On May 12, 2024, at 2:20 AM, a fire broke out in a laboratory at a university in Shanghai. The fire damaged some equipment, affecting an area of about 1m². Fortunately, there were no injuries reported.

Cause of the Incident

1. The laboratory's "oil bath" (a heated constant temperature magnetic stirrer) operated for an extended period, which leads to a failure in the temperature control system that caused the rising temperature and then a fire.
2. The experiment was left unattended overnight.

Important Reminders

1. Investigations revealed that the equipment involved was manufactured by certain companies. Laboratory personnel should be aware of the risks associated with this type of equipment and exercise caution in both purchasing and using it.

Collector type thermostatic
heating magnetic stirrer

2. Overnight experiments should be evaluated for potential risks and implement appropriate safety measures. Overnight experiments involving reactive, high-temperature, high-pressure, vacuum, and non-sealed Class 3B lasers, as well as running water require the completion of an "Overnight Experiment Application Form." This form must be signed by the PI and submitted to the laboratory safety officer and the safety director and then to HSE.
3. Safety Guidelines for Equipment Use:
 - Ensure the instrument is placed on a stable, clean, non-slip, dry, and fireproof surface, away from heat sources. Do not block the ventilation openings or heat dissipation areas.
 - Select suitable protective gear. When handling toxic or volatile substances, use closed containers and conduct all work within a suitable fume hood.
 - Before operating the equipment, set the speed control knob to the lowest setting.
 - Never allow the device to run dry! Exercise caution with high temperatures!
 - Regularly replace the heating medium in the “oil bath”. If the equipment will not be used for an extended period, drain any remaining water or oil, wipe it dry, unplug it, and store it in a cool, dry place.
 - If there are any unusual odors or sounds, immediately disconnect the power.
 - The layout of the laboratory should carefully consider the placement of equipment that carries higher risks, such as high-temperature and high-pressure devices. In labs with elevated operational risks, it may be beneficial to install monitoring systems and increase the frequency of training for staff working in those areas.

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